
pET-26b(+)

载体编号 Vector ID	V011007
抗性 Antibiotic Resistance	Kanamycin
长度 Length	5356 bp
载体类型 Vector Type	pET & Duet Vectors (Novagen)

产品信息

复制子 Replication Origin	ori
来源 Source/Author	Novagen (EMD Millipore)
拷贝数 Copy Number	High copy number
培养菌株 Growth Strain(s)	DH10B
培养温度 Growth Temperature	37℃
备注 Comment	pET-26b(+) is a prokaryotic expression vector with T7/lac promoter, pelB signal peptide for periplasmic secretion, C-terminal His-tag, Kan resistance, for secreted recombinant protein expression in E. coli.

质粒图谱



复溶方法

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add sterile water to dissolve the DNA: add 20 µl for 5 µg plasmid, and 100 µl for 100 µg plasmid.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 ℃.