
pCAGGS-ADRA1D(human)-EGFP-Strep-Tag II

载体编号 Vector ID	V025084
登录号 Accession ID	NM_000678.4
抗性 Antibiotic Resistance	Ampicillin
长度 Length	7268 bp
载体类型 Vector Type	Protein expression
宿主 Host	Mammalian cells

产品信息

复制子 Replication Origin	ori
筛选标记 Selection Marker	EGFP
启动子 Promoter	CAG
融合标签 Fusion Tag	EGFP-Strep-Tag II
表达方式 Expression Method	Transient

质粒图谱



复溶方法

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add sterile water to dissolve the DNA: add 20 µl for 5 µg plasmid, and 100 µl for 100 µg plasmid.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 ℃.