
pT3-EF1a-Kras(mouse)-3xFLAG-V5 tag

载体编号 Vector ID	V047648
登录号 Accession ID	NM_021284.7
抗性 Antibiotic Resistance	Ampicillin
长度 Length	6190 bp
载体类型 Vector Type	Protein expression
宿主 Host	Mammalian cells

产品信息

复制子 Replication Origin	ori
拷贝数 Copy Number	High
启动子 Promoter	EF1a
5' 引物 5' Primer	pEF-F
融合标签 Fusion Tag	3xFLAG-V5 tag
培养菌株 Growth Strain(s)	DH5a
培养温度 Growth Temperature	37℃
表达方式 Expression Method	Transient
备注 Comment	1 NT, 1 AA Mismatches

质粒图谱



复溶方法

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add sterile water to dissolve the DNA: add 20 µl for 5 µg plasmid, and 100 µl for 100 µg plasmid.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 ℃.